

### Mini Ball Implant 2.0 mm and 2.2 mm in 4 different lengths

- Auto-Advancing & Self-Tapping 2.0 mm and 2.2 mm diameter Mini Implant with 1.8 mm diameter Ball Attachment.
- Recommended for immediate loading procedures.
- Titanium Grade 5 Alloy Implant available in 10 mm, 11.5 mm, 13 mm and 15 mm lengths.
- Smooth and highly polished 2 mm collar.
- Surface prepared for osseous integration by Double Acid Etching and Grit Blasting.
- Implant provided in sterile packaging.
- Depending on bone density, drill 1/2 to 3/4 the length of implant.
- Purchase Mini Implant separately or with Mini Housing

#### BALL IMPLANT plus HOUSING

1 to 9 implants - \$65.00 each  
10 to 35 implants - \$62.50 each  
36 plus implants - \$58.00 each

#### BALL IMPLANT alone

1 to 9 implants - \$45.00 each  
10 to 35 implants - \$42.50 each  
36 plus implants - \$40.00 each

Ball Implant 2.0 mm Ø	
LENGTH	CAT. NO.
10 mm	EB2010
11.5 mm	EB2011
13 mm	EB2013
15 mm	EB2015

Ball Implant 2.2 mm Ø	
LENGTH	CAT. NO.
10 mm	EB2210
11.5 mm	EB2211
13 mm	EB2213
15 mm	EB2215



### Mini Square Implants 2.2 mm

- Auto-Advancing & Self-Tapping 2.2 mm diameter Implant with
- Preparable top recommended for immediate loading procedures.
- Titanium Grade 5 Alloy Implant available in 10 mm, 11.5 mm, 13 mm and 15 mm lengths.
- Smooth and highly polished 2 mm collar.
- Surface prepared for osseous integration by Double Acid Etching and Grit Blasting.
- Implant provided in sterile packaging.
- Depending on bone density, drill 1/2 to 3/4 the length of implant.

#### SQUARE IMPLANT

1 to 9 implants - \$55.00 each  
10 to 35 implants - \$50.00 each  
36 plus implants - \$45.00 each

Square Implant 2.2 mm Ø	
LENGTH	CAT. NO.
10 mm	EQ2210
11.5 mm	EQ2211
13 mm	EQ2213
15 mm	EQ2215

### Lance Drill 1.5 mm Ø

- Perforating the gingiva.

Cat. No. FL 1.5  
\$59.00 each



### Twist Drill 1.6 mm Ø

- Performing the osteotomy.

Cat. No. EIDL16  
\$59.00 each



### Stainless Steel Standard Housing 5.0 mm Ø with Silicone O-Ring

- Standard Housing with Silicone o-ring provides snug fit for 1.8 mm ball.

Cat. No. MH02A  
\$30.00 each



### Stainless Steel Mini Housing 3.8 mm Ø with Silicone O-Ring

- Mini Housing with Silicone o-ring provides snug fit for 1.8 mm ball.

Cat. No. MH03A  
\$30.00 each



### Silicone O-Rings for Standard Housing

Cat. No. MH02A-OR  
Pack of 10 O-Rings  
\$40.00 each



### Silicone O-Rings for Mini Housing

Cat. No. MH03A-OR  
Pack of 10 O-Rings  
\$40.00 each



### Implant Seating Tool Long

- Long shank used for seating mini implant in anterior arch.

Cat. No. EBTL  
\$90.00 each



### Implant Seating Tool Short

- Short shank used for seating mini implant in posterior arch.

Cat. No. EBTC  
\$90.00 each



### Implant Mount Driver

- For seating of implant using implant motor

Cat. No. IMDL  
\$75.00 each



### Thumb Driver

- For hand seating of implant.

Cat. No. ADL  
\$55.00 each



"Mini Implants are the most appreciated service that you can do for your patients. Stabilizing formerly intolerable dentures will truly change your patients' lives. Now they can eat, they can talk, and they can smile! With careful planning, the osteotomy can be done avoiding a flap resulting in a procedure almost without discomfort. The real satisfaction is when your patients leave your office with their stabilized dentures.  
Dr. Lloyd R. Taylor, President"



Each sealed individual implant is packaged in an outer vial.



Ask assistant to open the outer vial and drop the inner sterile vial with the implant onto your work area.



Open the inner sterile vial and remove Implant in its "ready to insert" placement holder.